

Getinge Atrium Chest Drain Allocation & Supply Update

Dear Valued GPO Partner:

We would like to provide you with an update regarding the allocation and supply for Atrium Chest Drain product lines.

Supply limitations related to COVID-19 and the global ABS plastic supply shortage continue to be managed, as they affect raw materials and components used in the manufacturing of the Atrium Chest Drains.

Based on this ongoing challenge, Getinge has stopped fulfilling Distributor orders in anticipation of taking orders directly from our mutual end-user customers. This action has been employed in an effort to more efficiently manage the allocation process and expedite shipping these orders.

We will continue to experience limited material availability through at least December, depending on when we can get increases in the supply of these raw materials and components. This has affected the Ocean, Oasis, Express, Pneumostat, and Express Mini 500 chest drain products. Intermittent availability will vary among all product codes. Once our supply improves, we will notify our distributor partners and revert back to the normal distribution process.

We are requesting that your member facilities place orders for full case quantities only, directly with Getinge and not their distributor for the following affected chest drain product codes:

Product Code	Product Description
16100	GETINGE PNEUMOSTAT CHEST DRN VALVE 10/CS
16400	GETINGE EXPRESS MINI 500 CHEST DRN 6/CS
2002-000	GETINGE OCEAN DRAIN, SINGLE W/AC,S 6/CS
2002-100	GETINGE OCEAN DRAIN, SINGLE W/O S 6/CS
2002-300	GETINGE OCEAN DRAIN, SINGLE W/S 6/CS
2002-400	GETINGE OCEAN DRAIN, SINGLE 6/CS
2012-320	GETINGE OCEAN DRAIN, PED W/AC,S 6/CS
2020-000	GETINGE OCEAN DRAIN, DUAL W/AC,S 6/CS
2020-300	GETINGE OCEAN DRAIN, DUAL W/S 6/CS
2050-000	GETINGE OCEAN DRAIN, BRU W/AC,S 6/CS
3600-100	GETINGE OASIS DRAIN, SINGLE W/AC 6/CS
3612-100	GETINGE OASIS DRAIN, PED W/AC 6/CS
3620-100	GETINGE OASIS DRAIN, DUAL AC 6/CS
3650-100	GETINGE OASIS DRAIN, BRU W/AC 6/CS
4000-100N	GETINGE EXPRESS DRAIN, SINGLE W/AC 6/CS
4050-100N	GETINGE EXPRESS DRAIN, BRU W/AC 6/CS

As an additional resource, attached is a cross-reference list of recommended product substitutes to help meet your members' needs for thoracic drainage products. Please note, other vendors may be experiencing supply disruptions as well.

Our commitment remains to provide you, our customers with reliable products and service. We apologize for any disruption this has caused.

Please contact Getinge Customer Service at (888) 943-8872 / (888-9GETUSA) or via email at vsPurchaseOrders@getinge.com to place an order or if you have additional questions.

Sincerely,



Charles Merchant PA
Sr. Director
Global Therapy Development
Acute Care Therapies

Getinge-to-Getinge Substitutes:

Chest Drains	Getinge Model Number Coding:
Ocean	2002 (Single), 2020 (Dual), 2012 (Pediatric), 2050 (ATS)
Oasis	3600 (Single), 3620 (Dual), 3612 (Pediatric), 3650 (ATS)
Express	4000 (Single), 4050 (ATS)

Chest Drains: If the backorder is for:	Recommended Getinge substitute:
Dual Collection Chest Drain	Two Single Collection Drains of the same model (i.e. Ocean, Oasis)
2020 3620	(2) 2002 (2) 3600
ATS Blood Recovery Chest Drain	Single Collection Drain of the same model (i.e. Ocean, Oasis, or Express)
2050 3650 4050	If no ATS being performed, then 2002 If no ATS being performed, then 3600 If no ATS being performed, then 4000
Single Collection Chest Drain	Recommend staying within the same model (i.e. Ocean, Oasis, or Express) and go to a Dual Collection Drain or ATS.
2002 3600	2020 or 2050 3620 or 3650
Pediatric Collection Chest Drain	Single Collection Drain of the same model (i.e. Ocean or Oasis) + Accessories
2012 3612	2002 + 19915 (16-40Fr) + 19918 (8, 10, 12Fr) 3600 + 19915 (16-40Fr) + 19918 (8, 10, 12Fr)
Pneumostat Chest Drain Valve	Cook Chest Drain Valve; Model # G36370
Express Mini 500 Chest Drain	Teleflex Mini Sahara; Model # S-0500

Getinge-to-Non Getinge Substitutes:

Single Chest Drain Series	Other recommended substitutes:
Ocean 2002	Teleflex: <i>Pleur-evac</i> A7000-08L, A7000-08LF, A8000-08L or A8000-08LF Cardinal: <i>Aqua Seal</i> 8888571315 or 8888571299 <i>Thora Seal 3</i> 8884713176 or 8884713308
Oasis 3600	Teleflex: <i>Pleur-evac</i> A6000-08L or A6000-08LF <i>Sahara</i> S1100-08L or S1100-08LF Cardinal: <i>Sentinel Seal</i> 8888571489 or 8888571562 <i>Altitude</i> 8888571370
Express 4000	Teleflex: <i>Sahara</i> S1100FS, S1100-08L, S1100-08LF or S1100-08FS
Pediatric Chest Drain Series	Other recommended substitutes:
Ocean 2012	Teleflex: <i>Pleur-evac</i> A8020-08L, A8020-08LF
Oasis 3612	Teleflex: <i>Pleur-evac</i> A6020-08L, A6020-08LF <i>Sahara</i> S1130-08LF
Dual Chest Drain Series	Other recommended substitutes:
Ocean 2020	Teleflex: <i>Pleur-evac</i> A7002-08L, A7002-08LF, A8002-08L or A8002-08LF
Oasis 3620	Teleflex: <i>Pleur-evac</i> A6002-08L, A6002-08LF <i>Sahara</i> S1102-08LF Cardinal: <i>Sentinel Seal</i> 8888571513
Blood Recovery Unit (ATS) Chest Drain Series	Other recommended substitutes:
Ocean 2050	Teleflex: <i>Pleur-evac</i> A9350-08L, A9350-08LF
Oasis 3650	Teleflex: <i>Pleur-evac</i> A9250-08L or A9250-08LF <i>Sahara</i> S1150-08L or S1150-08LF
Express 4050	Teleflex: <i>Sahara</i> S1150-08L or S1150-08LF