



PROFECT[®] HP

Hydrogen Peroxide Disinfectant



Use Profect HP daily as part of a simple and effective cleaning and disinfection program for your entire facility on hard, non-porous surfaces. Featuring patented hydrogen peroxide technology, Profect HP kills bacteria and viruses in 60 seconds. One minute contact times ensure efficacy and compliance for your most critical disinfection needs. Ideal for daily use on high-touch surfaces, Profect HP is available in a convenient ready-to-use quart. Effective, easy and safe, use Profect HP for all your one-step disinfecting, deodorizing, and cleaning needs.

EPA Registration Number: 5741-33

PROFECT® HP

FEATURES

- Kills bacteria, enveloped and non-enveloped viruses, and fungi
- Features patented hydrogen peroxide technology
- Ready-to-use formula ensures proper ppm required for disinfection
- Kills SARS-CoV-2 in 30 seconds

BENEFITS


- One minute contact times ensure efficacy and compliance
- Offers both safety and broad-spectrum kill in one product
- Meets the surface disinfection recommendations of OSHA bloodborne pathogen standard
- Disinfects high-touch surfaces in 1 minute

SPECIFICATION DATA

Dilution	Ready-to-use
pH	1.8 – 2.2
Color	Colorless
Fragrance	Mild
Specific Gravity @ 24°C/75°F	1.014
Flash Point	> 94°C/201°F ASTM D56
Stability: Shelf @ 24°C/75°F	1 year minimum
Stability: Freeze/Thaw	Will withstand 3 cycles
EPA Registration Number	5741-33

ADDITIONAL RESOURCES

Scan the QR code to access this product's safety information and documentation.



EFFICACY DATA

FOR USE AS A ONE-STEP DISINFECTANT, DEODORIZER, AND CLEANER ON HARD, INANIMATE NONPOROUS ENVIRONMENTAL SURFACES

When used as directed, this product is effective against *Pseudomonas aeruginosa*, *Salmonella enterica*, *Staphylococcus aureus*, *Acinetobacter baumannii*, *Escherichia coli* 0157:H7, Methicillin Resistant *Staphylococcus aureus* (MRSA), Multi-drug Resistant (MDR) *Enterococcus faecium*.

TUBERCULOCIDAL ACTIVITY

Profect HP exhibits disinfectant efficacy against *Mycobacterium bovis* (TB) when the treated surface is allowed to remain wet for 3 minutes at 68°F (20°C).

FUNGICIDAL ACTIVITY

Profect HP is an effective disinfectant against the pathogenic fungus, *Trichophyton interdigitale* (Athlete's Foot Fungus) (a cause of ringworm) on visibly clean, hard, non-porous surfaces (such as those found in locker rooms, dressing rooms, shower and bath areas, and exercise areas) for a 1-minute contact time at room temperature.

*VIRUCIDAL ACTIVITY ON HARD, INANIMATE NONPOROUS ENVIRONMENTAL SURFACES

When used as directed, Profect HP is effective against Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), Human Immunodeficiency Virus Type 1 (HIV-1), Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Human Coronavirus, Human Rotavirus, Influenza A Virus, Influenza B Virus, and Norovirus when treated surface is allowed to remain wet for 1 minute.

When used as directed, Profect HP is effective against Severe Acute Respiratory Syndrome-related Coronavirus 2 (SARS-COV-2) in 30 seconds!

ORDERING INFORMATION



#100803



SAFETY: For institutional and industrial use only. Be sure to read all directions, precautionary and first aid statements on product labels before using this or any other Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor and at www.spartanchemical.com. Workplace labels are available in English, Spanish, and French.

GUARANTEE: Spartan's modern manufacturing and laboratory control ensures uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within 12 months of the date of manufacture.

www.spartanchemical.com • 1-800-537-8990 • 1110 Spartan Drive, Maumee, OH 43537

© 2021 Spartan Chemical Company, Inc. All Rights Reserved. Profect is a registered trademark of Spartan Chemical Company, Inc.



10/21 L1008